


## Fluorescein II CEP (6-FAM II, BA 0253)

### Product Information



*Installs a 6-carboxyfluorescein internally or at the 5'-terminus of an oligonucleotide using a DMT-bearing phosphoramidite, allowing DMT-on cartridge purification.*

For the installation of fluorescein at the 5'-terminus of an oligonucleotide, the phosphoramidite "6-FAM" (5'-Fluorescein CEP, BA 0054), which does not bear a DMT group, is a popular choice. However, the lack of a trityl group precludes multiple additions or assaying the coupling step. We therefore offer Fluorescein II CEP ("6-FAM II", BA 0253), which features the same tether length as 6-FAM, but includes a DMT group.<sup>1</sup>

**Coupling:** Couple using normal instrument protocols except the coupling time should be extended to 15 min. Typical coupling yields are ca. 95%. Trityl-on mode is recommended (see below).

**Cleavage and nucleobase deprotection:** Use standard techniques. For the highest yields, prepare the oligo DMT-on and remove the trityl group after cleavage and deprotection. The DMT may also be used to facilitate cartridge purification with on-column detritylation, e.g. with Fluoro-Pak columns (see below).

**Purification:** Standard methods are applicable.

Fluoro-Pak™ columns may be used for cartridge purification with on-column detritylation if the final DMT group is left on. Use the following protocol:

#### 1. Materials needed:

- DMT-on oligonucleotide with 5'-fluorescein II installed. Ammonia may be present.
- 1 Fluoro-Pak™ or Fluoro-Pak™ II column (FP 7210 or FP 7220)
- 1 Luer adaptor (available with columns)
- 1 3 mL PE/PP syringe with male Luer tip
- 1 20 mL vial for catching waste
- 2 mL acetonitrile (MeCN)
- 3 mL 0.1 M aqueous TEAA (triethylammonium acetate)
- 2-5 mL Loading Buffer (Berry & Associates # LB-7100)
- 1 mL 3% aqueous TFA (trifluoroacetic acid)
- 1 mL 5% MeCN in 0.1 M aqueous TEAA
- 1 mL 50% MeCN in water

## 2. Sample preparation

Without removing the ammonia used in the deblocking step, dilute the crude deprotected oligonucleotide with an equal volume of Loading Buffer. The final sample volume should be 2-6 mL.

## 3. Conditioning the column

Pass the following through the Fluoro-Pak™ column to waste at a flow rate of 2 seconds per drop:

1. 2 mL of MeCN
2. 2 mL of 0.1 M aqueous TEAA
3. 2 mL of Loading Buffer

The first 2-3 drops of MeCN may be discolored.

## 4. Loading the DMT-on 6-FAM II-labeled oligonucleotide

Using the sample as prepared above, pass it through the pre-conditioned column at a flow rate of **5 seconds per drop**, conditions that allow one-pass loading. Most of the failure sequences pass through the column during loading.

## 5. Eluting remaining failure sequences

Pass the following through the Fluoro-Pak™ column to waste at a flow rate of **2 seconds per drop**:

1. 1 mL 5% MeCN in 0.1 M aqueous TEAA
2. 1 mL of water

The acetonitrile/TEAA solution selectively washes the remaining failure sequences from the resin while the fluorescein-labeled oligonucleotide is retained. The water step washes the buffer out of the resin prior to the acid-catalyzed detritylation in the next step.

## 6. On-column detritylation

Pass the following through the Fluoro-Pak™ column to waste at a flow rate of **2 seconds per drop**:

1. 2 mL 3% aqueous TFA. Column may become orange.
2. 1 mL of 0.1 M aqueous TEAA.
3. 1 mL of water.

## 7. Elution of the final detritylated oligonucleotide

Pass the following through the Fluoro-Pak™ column at **2 seconds per drop**, collecting the eluate in an appropriate sample tube:

1 mL of 50% acetonitrile in water. The oligo should elute in the first 15-20 drops.

Determine the optical density units at 494 nm to quantify the final oligonucleotide. Use as-is or lyophilize for short-term storage. If the oligonucleotide needs to be stored, add 100 mL of 10x TE buffer (100 mM Tris, 10 mM EDTA, pH 8) and store at 4 °C for a ready-use solution, or freeze at -20 °C for longer periods.

---

## References

1. 5'-Fluorescein II CEP incorporates 6-carboxyfluorescein. The 5-carboxyfluorescein isomer of this product is also known and has been used with cartridge purification using OPC columns. See: Theisen, P.; McCollum, C.; Upadhy, K.; Jacobson, K.; Vu, H.; Andrus, A. *Tetrahedron Lett.* **1992**, *33*, 5033-5036.