


BlackBerry[®] Quencher 650 CPG (3'-BBQ-650[®] CEP)
Product No. BL 2010
Product Information



Standard phosphoramidites and synthesis protocols may be employed with 3'-BBQ-650[®] CPG. Nucleobase deprotection may be performed using standard protocols with ammonium hydroxide or AMA without degradation of the quencher, e.g., 12 h, 55 °C. See Summary below for comparison to BL 2020 and BL 2030.

The lipophilicity of the BBQ-650[®] moiety may require the use of relatively high concentrations of the organic mobile phase in RP-HPLC purifications, especially with shorter oligonucleotides.

For quantification, the following extinction coefficients may be useful, which were determined using a simple BBQ-650[®] chromophore (i.e., no oligonucleotide): At 598 nm in methanol, $\epsilon = 40,667 \text{ M}^{-1}\text{cm}^{-1}$; at 260 nm in methanol, $\epsilon = 15,077 \text{ M}^{-1}\text{cm}^{-1}$.

3'-BBQ-650[®] CPG Product Comparison Summary

Compound	Column Wash	Cleavage Conditions	Deprotection
3'-BBQ-650[®] CPG (BL 2010)	None	Ammonium Hydroxide or AMA; 2 hours	12 hours, 55 °C
3'-BBQ-650[®] CPG II (BL 2020)	10% Diethyl Amine for 5-10 minutes, then blow dry 10 min.	Ammonium Hydroxide; 5 min. or AMA; 2 hours	65 °C for 10 min.
3'-BBQ-650[®] CPG III (BL 2030)	None	AMA only; 20 min.	65 °C for 10 min.